






Which SenseAnywhere datalogger do I need?



LN2	-80°C freezer Dry Ice shipment	Freezer Fridge Ambient Incubator	Oven Autoclave	Oven	High temp oven	Sensor type needed
						20-20-31 or 20-20-43 + 01-01-30 or 20-20-41 + Magnetic Door Contact Set
						20-20-43 + 01-01-00
						20-20-43 + 01-01-03 20-20-43 + 01-01-04
						20-20-43 + 01-01-02 20-20-43 + any 4-wire Pt100
						20-20-25 20-20-43 + 01-01-20
						20-20-30 (In transport packaging or room monitoring – shorter range)
						20-20-24 (Temp + RH in transport packaging or room monitoring – shorter range)

- Notes:**
- Besides the data loggers you will at least need one AccessPoint (infrastructure), depending on the area which you would like to cover multiple AccessPoints might be needed.
 - If you already have a running installation and you do not want to extend the wireless coverage you just add data loggers to your system.
 - If you need initial calibration and regular re-calibration we strongly recommend using an AssetTag for each device that needs to be calibrated. The system will keep track of the logger swapping without any downtime or administration. Keeping all history linked to the asset that is being monitored without paperwork.

EAN Codes

Product	Model	EAN EUR+Africa	EAN Americas	EAN US Pharma Grade
AiroSensor	20-20-24	0630125977842	0630125977941	0096647495200
AiroSensor ER	20-20-25	0630125977859	0630125977958	0096647495217
AiroSensor T	20-20-30	0630125977422	0630125977644	0096647495224
AiroSensor T ER	20-20-31	0630125977439	0630125977651	0096647495231
AiroSensor X ER	20-20-41	0630125977552	0630125977668	0096647495248
AiroSensor X ER	20-20-43	0096647494807	0096647494814	0096647495255
Indoor AccessPoint EU	20-01-10	0630125977064	-	-
Indoor AccessPoint UK	20-01-10	0630125977071	-	-
Indoor AccessPoint US	20-01-10	-	0630125977699	-
Mobile AccessPoint GPS 12/24	20-01-10	0630125977118	-	-
Mobile AccessPoint GPS 230	20-01-10	0630125977279	-	-
Mobile AccessPoint without GPS 12/24	20-01-10	0096647494944	0096647495286	-
Mobile AccessPoint without GPS 230	20-01-10	0096647494951	-	-
Mobile AccessPoint without GPS US 120	20-01-10	-	0096647495293	-
Outdoor AccessPoint Passive EU	20-01-10	0630125977088	-	-
Outdoor AccessPoint Active Un	20-01-10	0096647495132	-	-
Outdoor AccessPoint Passive UK	20-01-10	0630125977019	-	-
Outdoor AccessPoint Passive US	20-01-10	-	0630125977705	-
Outdoor AccessPoint Active US	20-01-10	-	0096647495156	-
Pt100 Probe	01-01-00	0096647494821	0096647494821	-
Pt100 Interface	01-01-02	0096647494838	0096647494838	-
Hermetically Sealed Pt100 5M	01-01-03	0096647495095	0096647495095	-
Hermetically Sealed Pt100 10M	01-01-04	0096647495163	0096647495163	-
Temp + RH Module	01-01-20	0096647495040	0096647495040	-
Temp Only Module	01-01-30	0096647495101	0096647495101	-
Magnetic Door Contact Set	-	0096647495187	0096647495187	-

Product Pictures



20-20-24
20-20-30



20-20-25
20-20-31



20-20-41
20-20-43



Indoor
AccessPoint



Mobile
AccessPoint With GPS



Mobile
AccessPoint Without GPS



Outdoor AccessPoint
Active/Passive



Pt100
Probe



Pt100
Hermetically Sealed



Pt100
Interface



Temp Only
Temp + RH Module



Magnetic Door
Contact Set